


Product Name: Recombinant Human CNKSR3 Protein
Catalog #: PEH7762



Summary

Name	CNKR3_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q6P9H4
Host	E.coli
Species	Human
Molecular Mass	37 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	QVVDWTRGLD DCLQQYVHKF EREKINGEQL LQISHQDLEE LGVTRIGHQE LVLEAVDLLC ALNYGLETDN MKNLVLKLRA SSHNLQNYIS SRRKSPAYDG NTSRKAPNEF LTSVVELIGA AKALLAWLDR APFTGITDFS VTKNKIIQLC LDLTTTVQKD CFVAEMEDKV LTVVKVLNGI CDKTIRSTTD PVMSQCACLE EVHLPNIKPG EGLGMYIKST YDGLHVITGT TENSADRSQ KIHAGDEVIQ VNQQTVVGWQ LKNLVKKLRE NPTGVVLLK KRPTGSFNFT PAPLKNLRWK PPLVQTSPPP

image



Background

Alternative Names	Connector enhancer of kinase suppressor of ras 3 (Connector enhancer of KSR 3) (CNK homolog protein 3) (CNK3) (CNKSR family member 3) (Maguin-
--------------------------	--

Product Name: Recombinant Human CNKSR3 Protein
Catalog #: PEH7762



like protein)

Note

For Research Use Only , Not for Diagnostic Use