


Product Name: Recombinant Human SKIC3 Protein
Catalog #: PEH7763



Summary

Name	TTC37_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q6PGP7
Host	E.coli
Species	Human
Molecular Mass	28 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	KSARDAIRNK EYKEALKHCK TVLKQEKNNY NAWVFIGVAA AELEQPDQAA SAYKKAEELE PDQLLAWQGL ANLYEKYNHI NAKDDLPGVY QKLLDLYESV DKQKWCDVCK KLVLDYYQEK KHLEVARTWH KLIKTRQEQG AENEELHQLW RKLTQFLAES TEDQNNETQQ LLFTAFENAL GLSDKIPSED HQVLYRHFIQ SLSKFPHESA RLKKACEGMI

image



Background

Alternative Names Tetratricopeptide repeat protein 37 (TPR repeat protein 37) (SKI3 homolog) (Ski3) (Tricho-hepatic-enteric syndrome protein) (Thespin)

Note

Product Name: Recombinant Human SKIC3 Protein
Catalog #: PEH7763



For Research Use Only , Not for Diagnostic Use