


Product Name: Recombinant Human PACS1 Protein
Catalog #: PEH7765



Summary

Name	PACS1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-MBP & C-His
Accession #	Q6VY07
Host	E.coli
Species	Human
Molecular Mass	73 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	EGIKSKLSDR SPDIDNYSEE EEESFSSEQE GSDDPLHGQD LFYEDEDLRK VKKTRRKLTS TSAITRQPNI KQKFVALLKR FKVSDEVGFG LEHVSREQIR EVEEDLDELY DSLEMYNPSD SGPEMEETES ILSTPKPKLK PFFEGMSQSS SQTEIGSLNS KGSLGKDTTS PMELAALEKI KSTWIKNQDD SLTETDTLEI TDQDMFGDAS TSLVVPEKVK TPMKSSKTDL QGSASPSKVE GVHTPRQKRS TPLKERQLSK PLSERTNSSD

image



Background

Alternative Names	Phosphofurin acidic cluster sorting protein 1 (PACS-1)
--------------------------	--

Product Name: Recombinant Human PACS1 Protein
Catalog #: PEH7765



Note

For Research Use Only , Not for Diagnostic Use