


Product Name: Recombinant Human MPRIP Protein
Catalog #: PEH7766



Summary

Name	MPRIP_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q6WCQ1
Host	E.coli
Species	Human
Molecular Mass	44 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	KFQANIFNKS KCQNCFKPRE SHLLNDEDLT QAKPIYGGWL LLAPDGTDFD NPVHRSRKWQ RRFILYEHG LLRYALDEMP TTLPQGTINM NQCTDVVDGE GRTGQKFSLC ILTPEKEHFI RAETKEIVSG WLEMLMVYPR TNKQNQKKKR KVEPPTQEP GPAKVAVTSS SSSSSSSSI PSAEKVPTTK STLWQEEMRT KDQPDGSSL PAQSPSQSQP PAASSLREPG LESKEESAM SSDRMDCGRK VRVESGYFSL EGTKQDLKAE EQQLPPPLSP PSPSTPNHRR SQVIEKFEAL DIEKAEHMET NAVGPSPSSD TRQGRSEKRA FPRKRDFTE APPAPLPDAS ASPLSPHRR KSLDRRSTEP

image



Background

Alternative Names Myosin phosphatase Rho-interacting protein (M-RIP) (Rho-interacting

Product Name: Recombinant Human MPRIP Protein
Catalog #: PEH7766



protein 3) (RIP3) (p116Rip)

Note

For Research Use Only , Not for Diagnostic Use