


Product Name: Recombinant Human BZW1 Protein
Catalog #: PEH7771



Summary

Name	BZW1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q7L1Q6
Host	E.coli
Species	Human
Molecular Mass	49 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	MNNQKQKPT LSGQRFKTRK RDEKERFDPT QFQDCIIQGL TETGTDLEAV AKFLDASGAK LDYRRYAETL FDILVAGGML APGGTLADDM MRTDVCVFAA QEDLETMQAF AQVFNKLIRR YKYLEKGFED EVKKLLLFLK GFSESERNKL AMLTGVLLAN GTLNASILNS LYNENLVKEG VSAFAVKLF KSWINEKDIN AVAASLRKVS MDNRLMELFP ANKQSVEHFT KYFTEAGLKE LSEYVRNQQT IGARKELQKE LQEQMSRGDP FKDIILYVKE EMKKNNIPEP VVIGIVWSSV MSTVEWNKKE ELVAEQAIKH LKQYSPLLA FTTQGOSELT LLLKIQEYCY DNIHFMKAFQ KIVVLFYKAE VLSEEPILKW YKDAHVAKGK SVFLEQMKKF VEWLKNAEEE SESEAEEGD

image



Background

Product Name: Recombinant Human BZW1 Protein
Catalog #: PEH7771



Alternative Names

Basic leucine zipper and W2 domain-containing protein 1 (Protein Orf)

Note

For Research Use Only , Not for Diagnostic Use