



Summary

Name	TAOK1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q7L7X3
Host	E.coli
Species	Human
Molecular Mass	54 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs.Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	LDMMEGDHTV MSNSSVIHLK PEEENYREEG DPRTRASDPQ SPPQVSRHKS HYRNREHFAT IRTASLVTRQ MQEHEQDSEL REQMSGYKRM RRQHQQQLMT LENKLKAEMD EHRLRLDKDL ETQRNNFAAE MEKLIKKHQA AMEKEAKVMS NEEKKFQQHI QAQQKELNS FLESQKREYK LRKEQLKEEL NENQSTPKKE KQEWLSKQKE NIQHFQAEED ANLLRRQRQY LELECRRFKR RMLLGRHNLE QDLVREELNK RQTQKDLEHA MLLRQHESMQ ELEFRHLNTI QKMRCELIRL QHQTELTNQL EYNKRREREL RRKHVMEVRQ QPKSLKSKEL QIKKQFQDTC KIQTRQYKAL RNHLLLETPK SEHKAVLKRL KEEQTRKLAI LAEQYDHSIN EMLSTQALRL DEAQEAECQV

image



Background

Product Name: Recombinant Human TAOK1 Protein
Catalog #: PEH7773



Alternative Names

Serine/threonine-protein kinase TAO1 (EC 2.7.11.1) (Kinase from chicken homolog B) (hKFC-B) (MARK Kinase) (MARKK) (Prostate-derived sterile 20-like kinase 2) (PSK-2) (PSK2) (Prostate-derived STE20-like kinase 2) (Thousand and one amino acid protein kinase 1) (TAOK1) (hTAOK1)

Note

For Research Use Only , Not for Diagnostic Use