


Product Name: Recombinant Human RASGRP2 Protein
Catalog #: PEH7774



Summary

Name	GRP2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q7LDG7
Host	E.coli
Species	Human
Molecular Mass	13 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	FQESNSLRPV ACRHCKALIL GIYKQGLKCR ACGVNCHKQC KDRLSVECRR RAQSVSLEGS APSPSPMHS HHRAFSFLP RRGRRSRPP

image



Background

Alternative Names	RAS guanyl-releasing protein 2 (Calcium and DAG-regulated guanine nucleotide exchange factor I) (Ca/DAG-GEFI) (Cdc25-like protein) (hCDC25L) (F25B3.3 kinase-like protein)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use