


Product Name: Recombinant Human ZC3HAV1 Protein
Catalog #: PEH7775



Summary

Name	ZCCHV_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q7Z2W4
Host	E.coli
Species	Human
Molecular Mass	21 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	SKHMQKNPPG PRAPSSHRRN MAYRARSKSR DRFFQGSQEF LASASASAER SCTPSPDQIS HRASLEDAPV DDLTRKFTYL GSQDRARPPS GSSKATDLGG TSQAGTSQRF LENGSEQEDLL HGPNPGSTYLA SNSTSAPNWK SLTSWTNDQG ARRKTVFSPT LPAARSSLGS LQTPEAVTTR

image



Background

Alternative Names	Zinc finger CCCH-type antiviral protein 1 (ADP-ribosyltransferase diphtheria toxin-like 13) (ARTD13) (Inactive Poly [ADP-ribose] polymerase 13) (PARP13) (Zinc finger CCCH domain-containing protein 2) (Zinc finger antiviral protein) (ZAP)
--------------------------	---

Product Name: Recombinant Human ZC3HAV1 Protein
Catalog #: PEH7775



Note

For Research Use Only , Not for Diagnostic Use