


Product Name: Recombinant Human HUWE1 Protein
Catalog #: PEH7777



Summary

Name	HUWE1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q7Z6Z7
Host	E.coli
Species	Human
Molecular Mass	46 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	FRQELERLDE GLRKEDMAVH VRRDHVFEDS YRELHRKSPE EMKNRLYIVF EGEEGQDAGG LLREWYMIIS REMFNPMYAL FRTSPGDRVT YTINPSSHCHN PNHLSYFKFV GRIVAKAVYD NRLLCYFTR SFYKHILGKS VRYTDMESD YHFYQGLVYL LENDVSTLGY DLTSTEVQE FGVCEVRDLK PNGANILVTE ENKKEYVHLV CQMRMTGAIR KQLAAFLEGF YEIIPKRLIS IFTEQELELL ISGLPTIDID DLKSNTEYHK YQSNSIQIQW FWRALRSFDQ ADRAKFLQFV TGTSKVPLQG FAALEGMNGI QKFQIHRDDR STDRLPSAHT CFNQLDLPAY ESFEKLRHML LLAIQECSEG

image



Background

Alternative Names	E3 ubiquitin-protein ligase HUWE1 (EC 2.3.2.26) (ARF-binding protein 1)
--------------------------	---

Product Name: Recombinant Human HUWE1 Protein
Catalog #: PEH7777



(ARF-BP1) (HECT, UBA and WWE domain-containing protein 1) (HECT-type E3 ubiquitin transferase HUWE1) (Homologous to E6AP carboxyl terminus homologous protein 9) (HectH9) (Large structure of URB1) (LASU1) (Mcl-1 ubiquitin ligase E3) (Mule) (Upstream regulatory element-binding protein 1) (URE-B1) (URE-binding protein 1)

Note

For Research Use Only , Not for Diagnostic Use