


Product Name: Recombinant Human ACOT1 Protein
Catalog #: PEH7781



Summary

Name	ACOT1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q86TX2
Host	E.coli
Species	Human
Molecular Mass	47 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	MAATLILEPA GRCCWDEPVR IAVRGLAPEQ PVTLRASLRD EKGALFQAHA RYRADTLGEL DLERAPALGG SFAGLEPMGL LWALEPEKPL VRLVKRDVRT PLAVELEVLD GHDPDPGRLL CRVRHERYFL PPGVRREPVR AGRVRGTLFL PPEPGPFPGI VDMFGTGGGL LEYRASLLAG KGFAVMALAY YNYEDLPKTM ETLHLEYFEE AVNYLLSHPE VKGPGVGLLG ISKGGECLCS MASFLKGITA AVVINGSVAN VGGTLRYKGE TLPPVGVNRRN RIKVTKDGYA DIVDVLNSPL EGPDQKSFIP VERAESTFLF LVGQDDHNWK SEFYANEACK RLQAHGRRKP QIICYPETGH YIEPPYFPLC RASLHALVGS PIIWGGEPRH HAMAQVDAWK QLQTFHKKHL GGHEGTIPSK V

image



Background

Product Name: Recombinant Human ACOT1 Protein
Catalog #: PEH7781



Alternative Names

Acyl-coenzyme A thioesterase 1 (Acyl-CoA thioesterase 1) (EC 3.1.2.-) (CTE-I)
(CTE-Ib) (Inducible cytosolic acyl-coenzyme A thioester hydrolase) (Long
chain acyl-CoA thioester hydrolase) (Long chain acyl-CoA hydrolase)
(Palmitoyl-coenzyme A thioesterase) (EC 3.1.2.2)

Note

For Research Use Only , Not for Diagnostic Use