



Product Name: Recombinant Human PPARGC1B Protein
Catalog #: PEH7790



Summary

Name	PRGC2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q86YN6
Host	E.coli
Species	Human
Molecular Mass	10 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	QAGLRSKSQR PCVKADSTQD KKAPMMQSQS RSCTELHKHL TSAQCCLQDR GLQPPCLQSP RLPKEDKEP

image



Background

Alternative Names	Peroxisome proliferator-activated receptor gamma coactivator 1-beta (PGC-1-beta) (PPAR-gamma coactivator 1-beta) (PPARGC-1-beta) (PGC-1-related estrogen receptor alpha coactivator)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use