


Product Name: Recombinant Human PHYKPL Protein
Catalog #: PEH7792



Summary

Name	AT2L2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q8IUZ5
Host	E.coli
Species	Human
Molecular Mass	51 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	MAADQRPKAD TLALRQLIS SSCRLFFPED PVKIVRAQQQ YMYDEQGAIEY IDCISNVAHV GHCHPLVVQA AHEQNQVLNT NSRYLHDNIV DYAQRLSETL PEQLCVFYFL NSGSEANDLA LRLARHYTGH QDVVLDHAY HGHLSLIDI SPYKFRNLDG QKEWVHVAPL PDTYRGPYRE DHPNPAMAYA NEVKRVVSSA QEKGRKIAAF FAESLPSVGG QIIPPAGYFS QVAEHIRKAG GVFVADEIQV GFGRVGKHFV AFQLQGKDFV PDIVTMGKSI GNGHPVACVA ATQPVARAFE ATGVEYFNTF GGSPVSCAVG LAVLNVLEKE QLQDHATSVG SFLMQLLGQQ KIKHPIVGDV RGVGLFIGVD LIKDEATRTP ATEEAAYLVS RLKENYVLLS TDGPGRNLIK FKPPMCFSLD NARQVVAKLD AILTDMEEKV RSCETLRLQP

image



Background

Product Name: Recombinant Human PHYKPL Protein
Catalog #: PEH7792



Alternative Names

5-phosphohydroxy-L-lysine phospho-lyase (EC 4.2.3.134) (Alanine--glyoxylate aminotransferase 2-like 2)

Note

For Research Use Only , Not for Diagnostic Use