


Summary

| | |
|------------------------------|--|
| Name | ANKY2_HUMAN |
| Purity | Greater than 85% as determined by reducing SDS-PAGE |
| Endotoxin level | Please contact with the lab for this information |
| Tag | C-His |
| Accession # | Q8IV38 |
| Host | E.coli |
| Species | Human |
| Molecular Mass | 50 kDa |
| Buffer | Lyophilized from Phosphate buffered saline |
| Shipping | The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below. |
| Stability&Storage | Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles. |
| Sequence | MVHIKKGELT QEEKELLEVI GKGTVQEAGT LLSSKNVRVN CLDENGMTPL MHAAYK GKLD MCKLLLRHGA DVNCHQHEHG YTALMFAALS GNKDITWVML EAGAETDVVN SVGR TAAQMA AFV GQHDCVT IINFFPRER LDYYTKPQGL DKEPKLPPKL AGPLHKIITT TNLHPVKIVM LVNENPLLTE EAALNKCYRV MDLICEKCMK QRDMNEVLAM KMHYISCIFQ KCINFLKDGE NKLDTLIKSL LKGRASDGFP VYQEKIRES IRKFPYCEAT LLQQLVRSIA PVEIGSDPTA FSVLTQAITG QVGFVDVEFC TTCGEGKASK RCSVCKMVIY CDQTCQKTHW FTHKKICKNL KDIYEKQQL EAAKEKQREEN HGKLDVNSNC VNEEQPEAEV GISQKDSNPE DSGEGKKESEAELEGLQ DAPAGPQVSE E |

image



Background

Product Name: Recombinant Human ANKMY2 Protein
Catalog #: PEH7793



Alternative Names

Ankyrin repeat and MYND domain-containing protein 2

Note

For Research Use Only , Not for Diagnostic Use