


Product Name: Recombinant Human SLC9A9 Protein
Catalog #: PEH7794



Summary

Name	SL9A9_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q8IVB4
Host	E.coli
Species	Human
Molecular Mass	19 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	VGVDLDENLK EDPSSQHQA NNLDKNMTKA ESARLFMRWY SFDHKYLPKPI LTHSGPPLTT TLPEWCGPIS RLLTSPQAYG EQLKEDDVEC IVNQDELAIN YQEQASSPCS PPARLGLDQK ASPQTPGKEN IYEGDLGLGG YELKLEQTLG

image



Background

Alternative Names	Sodium/hydrogen exchanger 9 (Na ⁺)/H ⁺ exchanger 9) (NHE-9) (Solute carrier family 9 member 9)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use