



Product Name: Recombinant Human KSR1 Protein
Catalog #: PEH7795



Summary

Name	KSR1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q8IVT5
Host	E.coli
Species	Human
Molecular Mass	13 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	HNKCTKEAPA CRISFLPLTR LRRTESVPSD INNPVDRAAE PHFGTLPKAL TKKEHPPAMN HLDSSSNPSS TTSSTPSSPA PFPTSSNPSS ATTPPNPSPG

image



Background

Alternative Names Kinase suppressor of Ras 1 (EC 2.7.11.1)

Note

For Research Use Only , Not for Diagnostic Use