


Product Name: Recombinant Human PARP9 Protein
Catalog #: PEH7796



Summary

Name	PARP9_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-SUMO & C-His
Accession #	Q8IXQ6
Host	E.coli
Species	Human
Molecular Mass	28 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	TREKRENGL EARSPLAINM GFNVEEMYEA HAWIQRILSL QNHIIENNH ILYLGRKEHD ILSQLQKTSS VSITEIISPG RTELEIEGAR ADLIEVVMNI EDMLCKVQEE MARKKERGLW

image



Background

Alternative Names Protein mono-ADP-ribosyltransferase PARP9 (EC 2.4.2.-) (ADP-ribosyltransferase diphtheria toxin-like 9) (ARTD9) (B aggressive lymphoma protein) (Poly [ADP-ribose] polymerase 9) (PARP-9)

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human PARP9 Protein
Catalog #: PEH7796

