


Product Name: Recombinant Human ALDH16A1 Protein
Catalog #: PEH7798



Summary

Name	A16A1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q8IZ83
Host	E.coli
Species	Human
Molecular Mass	48 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	PEHRNSVPCQ DPITGENLAS CLQAQAEDVA AAVEAARMAF KGWSAHPGVV RAQHLTRLAE VIQKHQRLW TLESVLTGRA VREVRDGDVQ LAQQLLHYHA IQASTQEEAL AGWEPMGVIG LILPPTFSFL EMMWRICPAL AVGCTVVALV PPASPAPLLL AQLAGELGPF PGILNVLSGP ASLVPILASQ PGIRKVAFCG APEEGRALRR SLAGECAELG LALGTESLLL LTDTADVDSA VEGVVDAAWS DRGPGGLRLL IQESVWDEAM RRLQERMGR LRSRGLDGAV DMGARGAAAC DLVQRFVREA QSQGAQVFQA GDVPSERPFY PPTLVSNLPP ASPCAQVEVP WPVVVASPFR TAKEALLVAN GTPRGGASV WSERLGQALE LYGGLQVGTV WINAHGLRDP SVPTGGCKES

image



Background

Product Name: Recombinant Human ALDH16A1 Protein
Catalog #: PEH7798



Alternative Names

Aldehyde dehydrogenase family 16 member A1

Note

For Research Use Only , Not for Diagnostic Use