


**Product Name: Recombinant Human ABI1 Protein**  
**Catalog #: PEH7799**



## Summary

<b>Name</b>	ABI1_HUMAN
<b>Purity</b>	Greater than 85% as determined by reducing SDS-PAGE
<b>Endotoxin level</b>	Please contact with the lab for this information
<b>Tag</b>	C-His
<b>Accession #</b>	Q8IZP0
<b>Host</b>	E.coli
<b>Species</b>	Human
<b>Molecular Mass</b>	13 kDa
<b>Buffer</b>	Lyophilized from Phosphate buffered saline
<b>Shipping</b>	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
<b>Stability&amp;Storage</b>	Store at $\leq -70^{\circ}\text{C}$ , stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
<b>Sequence</b>	THKIIAPANM ERPVRYIRKP IDYTVLDDVG HGVKWLKAKH GNNQPARTGT LSRTNPPTQK PPSPMSGRG TLGRNTPYKT LEPVKPPTVP NDYMTSPARL

## image



## Background

**Alternative Names** Abl interactor 1 (Abelson interactor 1) (Abi-1) (Abl-binding protein 4) (AbIBP4) (Eps8 SH3 domain-binding protein) (Eps8-binding protein) (Nap1-binding protein) (Nap1BP) (Spectrin SH3 domain-binding protein 1) (e3B1)

## Note

For Research Use Only , Not for Diagnostic Use