


**Product Name: Recombinant Human LRRC47 Protein**  
**Catalog #: PEH7800**



## Summary

<b>Name</b>	LRC47_HUMAN
<b>Purity</b>	Greater than 85% as determined by reducing SDS-PAGE
<b>Endotoxin level</b>	Please contact with the lab for this information
<b>Tag</b>	C-His
<b>Accession #</b>	Q8N1G4
<b>Host</b>	E.coli
<b>Species</b>	Human
<b>Molecular Mass</b>	39 kDa
<b>Buffer</b>	Phosphate buffered saline
<b>Shipping</b>	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
<b>Stability&amp;Storage</b>	Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
<b>Sequence</b>	RLEKMVSGCQ TRSILEYLRV GGRGGGKGGK RAEGSEKEES RRKRERKQR REGGDGEEQD VGDAGRLLLR VLHVSENPVP LTVRVSPEVR DVRPYIVGAV VRGMDLQPGN ALKRFLTSQT KLHEDLCEKR TAATLATHL RAVKGPLLYC ARPPQDLKIV PLGRKEAKAK ELVRQLQLEA EEQRKQKKRQ SVSGLHRYLH LLDGNENYPC LVDADGDVIS FPPITNSEKT KVKKTTSDLF LEVTSATSLQ ICKDVMDALI LKMAEMKKYT LENKEEGSLS DTEADAVSQG LPDPTTNPSA GKDGPSLLVV EQVRVVDLEG SLKVYPSKA DLATAPPHVT

## image



## Background

**Alternative Names** Leucine-rich repeat-containing protein 47

**Product Name: Recombinant Human LRRC47 Protein**  
**Catalog #: PEH7800**

---



**Note**

For Research Use Only , Not for Diagnostic Use