


Product Name: Recombinant Human KRT78 Protein
Catalog #: PEH7801



Summary

Name	K2C78_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q8N1N4
Host	E.coli
Species	Human
Molecular Mass	37 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	QFASFIDKVR FLEQQNKVLE TKWHLLQQQG LSGSQQGLEP VFEACLDQLR KQLEQLQGER GALDAELKAC RDQEEYKSK YEEAHRRAT LENDFVVLKK DVDGVFLSKM ELEGKLEALR EYLYFLKHLN EEELGQLQTQ ASDTSVVLMS DNNRYLDFSS IITEVRARYE EIARSSKAEA EALYQTKYQE LQVSAQLHGD RMQETKVQIS QLHQEIQLRQ SQTENLKKQN ASLQAAITDA EQRGELALKD AQAKVDELEA ALRMAKQNLA RLLCEYQELT STKLSLDVEI ATYRRLLGE

image



Background

Alternative Names	Keratin, type II cytoskeletal 78 (Cytokeratin-78) (CK-78) (Keratin-5b) (Keratin-78) (K78) (Type-II keratin Kb40)
--------------------------	--

Product Name: Recombinant Human KRT78 Protein
Catalog #: PEH7801



Note

For Research Use Only , Not for Diagnostic Use