



Product Name: Recombinant Human MT1M Protein
Catalog #: PEH7802



Summary

Name	MT1M_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-SUMO & C-His
Accession #	Q8N339
Host	E.coli
Species	Human
Molecular Mass	18 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	MDPNCSTTG VSACTGSCT CKECKTSCK KSCCSCPVG CAKCAHGVCV KGTLENCSCC A

image



Background

Alternative Names Metallothionein-1M (MT-1M) (Metallothionein-1M) (MT-1M)

Note

For Research Use Only , Not for Diagnostic Use