


Summary

| | |
|------------------------------|--|
| Name | RHG18_HUMAN |
| Purity | Greater than 85% as determined by reducing SDS-PAGE |
| Endotoxin level | Please contact with the lab for this information |
| Tag | C-His |
| Accession # | Q8N392 |
| Host | E.coli |
| Species | Human |
| Molecular Mass | 37 kDa |
| Buffer | Phosphate buffered saline |
| Shipping | The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below. |
| Stability&Storage | Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles. |
| Sequence | VLTAHPSGK DQTVGNHAK AGEEATSSRR YGQYTMNQES TTIKVMKPP FDRSISQDSL DELSMEDYWI ELENIKSSE NSQEDQEVVV VKEPDEGELE EEWLKEAGLS NLFGESAGDP QESIVFLSTL TRTQAAAVQK RVETVSQTLR KKNKQYQIPD VRDIFAQQRE SKETAPGGTE SQSLRTNENK YQGRDDEASN LVGEEKLIPP EETPAPETDI NLEVSFAEQA LNQKESSEK IQKSKGDDAT LPSFRLPKDK TGTRIGDLA PQDMKKVCHL ALIELTALYD VLGIELKQQK AVKIKTKDSG |

image



Background

| | |
|--------------------------|---|
| Alternative Names | Rho GTPase-activating protein 18 (MacGAP) (Rho-type GTPase-activating protein 18) |
|--------------------------|---|

Product Name: Recombinant Human ARHGAP18 Protein
Catalog #: PEH7804



Note

For Research Use Only , Not for Diagnostic Use