


Product Name: Recombinant Human SYNPO Protein
Catalog #: PEH7806



Summary

Name	SYNPO_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q8N3V7
Host	E.coli
Species	Human
Molecular Mass	29 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	GPTDLEEDG VSRSGDDSAC RVTQGTPQLP KALGIQPPSC SREEQGASQH DDRASQDWDV VKAGQMMTAS PSPGPGPRVA QKPALGRSTS LTKDLKEAK ARSQQIAAQL TTPSSNSRG VQLFNRRRQR VNEFTLESHG QRGQKPSQES LRVLPSSLPG HAPGLSLSST SLPEPGPPRH PSPQSPDRGV PGHSMEGYSE EASLLRHLEK VASEEEEEVPL VVYLKENAAL LTANGLHLSQ NREAQQSSPA

image



SDS

Background

Alternative Names	Synaptopodin
--------------------------	--------------

Note

Product Name: Recombinant Human SYNPO Protein
Catalog #: PEH7806



For Research Use Only , Not for Diagnostic Use