

Product Name: Recombinant Human OXR1 Protein
Catalog #: PEH7807



Summary

| | |
|------------------------------|---|
| Name | OXR1_HUMAN |
| Purity | Greater than 85% as determined by reducing SDS-PAGE |
| Endotoxin level | Please contact with the lab for this information |
| Tag | C-His |
| Accession # | Q8N573 |
| Host | E.coli |
| Species | Human |
| Molecular Mass | 38 kDa |
| Buffer | Phosphate buffered saline |
| Shipping | The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below. |
| Stability&Storage | Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles. |
| Sequence | KIKESLPIDI DQLSGRDFCH SKKMTGSNTE EIDSRIRDAG NDSASTAPRS TEELSSEDVF TEELSPIRE ELVSSDELRLQ DKSSGASSES VQTVNQAEVE SLTVKSESTG TPGHLRSDE HSTNEVGTLC HKTDLNNLEM AIKEDQIADN FQGISGPKED STSIKGNSDQ DSFLHENS LH QEESQKENMP CGETAEFKQK QSVNKGKQ GK EQNQDSQTEA EELRKLWKTH T MQQTKQQRE NIQQVSQKEA KHKITSADGH IESSALLKEK QRHRLHKFLC LRVGKPMRKT FVSQASATMQ QYAQRDKKHE YWFAVPQERT |

image



Background

Alternative Names Oxidation resistance protein 1

Product Name: Recombinant Human OXR1 Protein
Catalog #: PEH7807



Note

For Research Use Only , Not for Diagnostic Use