


Product Name: Recombinant Human MINDY1 Protein
Catalog #: PEH7808



Summary

Name	MINDY1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q8N5J2
Host	E.coli
Species	Human
Molecular Mass	52 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	MEYHQPEDPA PGKAGTAEAV IPENHEVLG PDEHPQDTDA RDADGEARER EPADQALLPS QCGDNLESPL PEASSAPPGP TLGTLPEVET IRACSMQPQL PQSPRTRQPE PDFYCVKWIP WKGEQTPIIT QSTNGPCPLL AIMNILFLQW KVKLPPQKEV ITSDELMAHL GNCLLSIKPQ EKSEGLQLNF QQNVDDAMTV LPKLATGLDV NVRFTGVSDF EYTPECSVFD LLGIPLYHGW LVDPQSPEAV RAVGKLSYNQ LVERIITCKH SSDTNLVTEG LIAEQFLETT AAQLTYHGLC ELTAAAKEGE LSVFFRNNHF STMTKHKSHL YLLVTDQGFL QEEQVVWESL HNVDGDSCFC DSDFHLSHSL GKGPGEAGGS GSPETQLQVD QDYLIASLQ QQQPRGPLGL TDLELAQQLQ QEEYQQQQA QPVRMRTRVL SLQGRGATSG RPAGERRQRP KHESDCILL

image



Background

Product Name: Recombinant Human MINDY1 Protein
Catalog #: PEH7808



Alternative Names

Ubiquitin carboxyl-terminal hydrolase MINDY-1 (EC 3.4.19.12)
(Deubiquitinating enzyme MINDY-1) (Protein FAM63A)

Note

For Research Use Only , Not for Diagnostic Use