


Product Name: Recombinant Human ARFGAP2 Protein
Catalog #: PEH7809



Summary

Name	ARFG2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q8N6H7
Host	E.coli
Species	Human
Molecular Mass	20 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	LGAQKVSSQS FSEIERQAQV AEKLREQQAA DAKKQAEESM VASMLAYQE LQIDRKKEEK KLQNLEGKKR EQAERLGMGL VSRSSVSHSV LSEMQUIEQE TPVSAKSSRS QLDFDDVGT FASGPPKYKD NPFSLGESFG SRWDTDAAWG MDRVEEKEPE

image



Background

Alternative Names	ADP-ribosylation factor GTPase-activating protein 2 (ARF GAP 2) (GTPase-activating protein ZNF289) (Zinc finger protein 289)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human ARFGAP2 Protein
Catalog #: PEH7809

