

Product Name: Recombinant Human BANP Protein
Catalog #: PEH7810



Summary

Name	BANP_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q8N9N5
Host	E.coli
Species	Human
Molecular Mass	57 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	MMSEHDLADV VQIAVEDLSP DHPVVLENHV VTDEDEPALK RQRLEINCQD PSIKTICLRL DSIEAKLQAL EATCKSLEEK LDLVTNKQHS PIQVPMVAGS PLGATQTCNK VRCVVPQTTV ILNDRQNAI VAKMEDPLSN RAPDSLENI SNAVPGRRQN TIVVKVPGQE DSHHEDGESG SEASDSVSSC GQAGSQSIGS NVTLITLNSE EDYPNGTWLG DENNPEMRVR CAIIPSDMLH ISTNCRTAEK MALTLLDYLF HREVQAVSNL SGQGKHGKKQ LDPLTIYGIR CHLFYKFGIT ESDWYRIKQS IDSKCRTAWR RKQRGQSLAV KSFSRTPNS SSYCPSEPM STPPASELP QPQPQPQALH YALANAQQVQ IHQIGEDGQV QVGHHLIAQV PQGEQVQITQ DSEGNLQIHH VGQDQQLLEA TRIPCLLAPS VFKASSGQVL QGAQLIAVAS SDPAAAGVDG SPLQGSDIQV QYVQLAPVSD HTAGAQTAEA LQPTLQPEMQ LEHGAIQIQ

image

SDS

Background

Alternative Names Protein BANP (BEN domain-containing protein 1) (Btg3-associated nuclear protein) (Scaffold/matrix-associated region-1-binding protein)

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human BANP Protein
Catalog #: PEH7810

