


Product Name: Recombinant Human TOR1AIP2 Protein
Catalog #: PEH7815



Summary

Name	TOIP2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q8NFQ8
Host	E.coli
Species	Human
Molecular Mass	24 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	EDSQKDLEND PSVNSQAQET TIIASNAEEA EILHSACGLS KDHQEVETEG PESADTGDKS ESPDEANVGK HPKDKTEDEN KQSFLDGGKG HHLPSENLGK EPLDPDPSHS PSDKVGRADA HLGSSSVALP KEASDGTGAS QEPPTTDSQE AQSPGHSSAG QEGEDTLRRR LLAPEAGSHP QQTQKLEEIK ENAQDTMRQI

image



Background

Alternative Names Torsin-1A-interacting protein 2 (Lumenal domain-like LAP1)

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human TOR1AIP2 Protein
Catalog #: PEH7815

