


Product Name: Recombinant Human CMAS Protein
Catalog #: PEH7816



Summary

Name	NEUA_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-MBP & C-His
Accession #	Q8NFW8
Host	E.coli
Species	Human
Molecular Mass	52 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	DSGAFQSVWV STDHDEIENV AKQFGAQVHR RSSEVSKDSS TSLDAIIIEFL NYHNEVDIVG NIQATSPCLH PTDLQKVAEM

image



Background

Alternative Names	N-acetylneuraminyl cytidyltransferase (EC 2.7.7.43) (CMP-N-acetylneuraminic acid synthase) (CMP-NeuNAc synthase)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use