


Summary

Name	CHD5_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q8TDI0
Host	E.coli
Species	Human
Molecular Mass	33 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	AISKLLDRNQ DATDDTELQN MNEYLSSFVKV AQYVVREEDG VEEVEREIIK QEENVDPDYW EKLLRHHEQ QQEDLARNLG KGKRIRKQVN YNDASQEDQE WQDELSDNQS EYSIGSEDED EDFEERPEGQ SGRRQSRRQL KSDRDKPLPP LLARVGGNIE VLGFNARQRK AFLNAIMRWG MPPQDAFNH WLVRDLRGKS EKEFRAYVSL FMRHLCEPGA DGAETFADGV PREGLSRQHV LTRIGVMSLV RKKVQEFHV

image



Background

Alternative Names	Chromodomain-helicase-DNA-binding protein 5 (CHD-5) (EC 3.6.4.12) (ATP-dependent helicase CHD5)
--------------------------	---

Product Name: Recombinant Human CHD5 Protein
Catalog #: PEH7820



Note

For Research Use Only , Not for Diagnostic Use