


**Product Name: Recombinant Human DLG5 Protein**  
**Catalog #: PEH7821**



## Summary

<b>Name</b>	DLG5_HUMAN
<b>Purity</b>	Greater than 85% as determined by reducing SDS-PAGE
<b>Endotoxin level</b>	Please contact with the lab for this information
<b>Tag</b>	C-His
<b>Accession #</b>	Q8TDM6
<b>Host</b>	E.coli
<b>Species</b>	Human
<b>Molecular Mass</b>	40 kDa
<b>Buffer</b>	Phosphate buffered saline
<b>Shipping</b>	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
<b>Stability&amp;Storage</b>	Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
<b>Sequence</b>	LPGDSFYIRA LYDRLADVEQ ELSFKDDIL YVDDTLPQGT FGSWMAWQLD ENAQKIQRGQ IPSKYVMDQE FSRRLSMSEV KDDNSATKTL SAAARRSFFR RKHKHKRSGS KDGKDLLALD AFSSDSIPLF EDSVSLAYQR VQKVDCTALR PVLILGPLLD VVKEMLVNEA PGKFCRCPLE VMKASQQAIE RGVKDCLFVD YKRRSGHFDV TTVASIKEIT EKNRHCLLDI APHAIERLHH MHIYPVIFI HYKSAKHIKE QRDPIYLRDK VTQRHSKEQF EAAQKLEQEQY SRYFTGVIQG GALSSICTQI LAMVNQEENK

## image



## Background

**Alternative Names** Disks large homolog 5 (Discs large protein P-dlg) (Placenta and prostate DLG)

**Product Name: Recombinant Human DLG5 Protein**  
**Catalog #: PEH7821**

---



**Note**

For Research Use Only , Not for Diagnostic Use