


Product Name: Recombinant Human DTD1 Protein
Catalog #: PEH7823



Summary

Name	DTD1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q8TEA8
Host	E.coli
Species	Human
Molecular Mass	24 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	MKAVVQRVTR ASVTVGGEQI SAIGRGICVL LGISLEDTQK ELEHMVRKIL NLRVFEDESG KHWSKSVMDK QYEILCVSQF TLQCVLKGNK PDFHLAMPTE QAEGFYNSFL EQLRKTYRPE LIKDGKFGAY MQVHIQNDGP VTIELESPAP GTATSDPKQL SKLEKQQQRK EKTRAKGPSE SSKERNTPRK EDRSASSGAE GDVSSEREP

image



Background

Alternative Names	D-aminoacyl-tRNA deacylase 1 (DTD) (EC 3.1.1.96) (DNA-unwinding element-binding protein B) (DUE-B) (Gly-tRNA(Ala) deacylase) (Histidyl-tRNA synthase-related)
--------------------------	---

Product Name: Recombinant Human DTD1 Protein
Catalog #: PEH7823



Note

For Research Use Only , Not for Diagnostic Use