


Product Name: Recombinant Human TTN Protein
Catalog #: PEH7829



Summary

Name	TITIN_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q8WZ42
Host	E.coli
Species	Human
Molecular Mass	33 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	GKDNWIRCNM KLVPELTYKV TGLEKGNKYL YRVSAENKAG VSDPSEILGP LTADDAFVEP TMDLSAFKDG LEVIVPNPIT ILVPSTGYPR PTATWCFGDK VLETGDRVKM KTLSAYAELV ISPERSDKG IYTLKLENRV KTISGEIDVN VIARPSAPKE LKFGDITKDS VHLTWEPPDD DGGSPLTGYV VEKREVSRT WTKVMDFVTD LEFTVPDLVQ GKEYLFKVCA RNKCGPGEPA YVDEPVMNST PATVPDPEN VKWRDRTANS IFLTWDPPKN

image



Background

Alternative Names Titin (EC 2.7.11.1) (Connectin) (Rhabdomyosarcoma antigen MU-RMS-40.14)

Product Name: Recombinant Human TTN Protein
Catalog #: PEH7829



Note

For Research Use Only , Not for Diagnostic Use