


Product Name: Recombinant Human LZIC Protein
Catalog #: PEH7830



Summary

Name	LZIC_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q8WZA0
Host	E.coli
Species	Human
Molecular Mass	22 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	MASRGKTETS KLKQNLEEQL DRLMQQLQDL EECREELDTD EYEETKKETL EQLSEFNDSL KKIMSGNMTL VDELSGMQLA IQAAISQAFK TPEVIRLFAK KQPGQLRTRL AEMDRDLMVG KLERDLYTQQ KVEILTALRK LGEKLTADDE AFLSANAGAI LSQFEKVSTD LGSGDKILAL ASFEVEKTKK

image



Background

Alternative Names	Protein LZIC (Leucine zipper and CTNNBIP1 domain-containing protein) (Leucine zipper and ICAT homologous domain-containing protein)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human LZIC Protein
Catalog #: PEH7830

