


Product Name: Recombinant Human HSD17B8 Protein
Catalog #: PEH7831



Summary

Name	DHB8_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-MBP & C-His
Accession #	Q92506
Host	E.coli
Species	Human
Molecular Mass	69 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	MASQLQNRLR SALALVTGAG SGIGRAVSVR LAGEGATVAA CDLDRAAAQE TVRLLGGPGS KEGPPRGNHA AFQADVSEAR AARCLLEQVQ ACFSRPPSVV VSCAGITQDE FLLHMSEDDW DKVIAVNLKG TFLVTQAAAQ ALVSNGCRGS IINISSIVGK VGNVGTNYA ASKAGVIGLT QTAARELGRH GIRCNSVLPG FIATPMTQKV PQKVVDKITE MIPMGHLGDP EDVADVVAFL ASED SGYITG TSVEVTGGLF M

image



Background

Alternative Names	Estradiol 17-beta-dehydrogenase 8 (EC 1.1.1.62) (17-beta-hydroxysteroid dehydrogenase 8) (17-beta-HSD 8) (3-ketoacyl-[acyl-carrier-protein] reductase alpha subunit) (KAR alpha subunit) (3-oxoacyl-[acyl-carrier-
--------------------------	--

Product Name: Recombinant Human HSD17B8 Protein
Catalog #: PEH7831



protein] reductase) (EC 1.1.1.-) (Protein Ke6) (Ke-6) (Really interesting new gene 2 protein) (Short chain dehydrogenase/reductase family 30C member 1) (Testosterone 17-beta-dehydrogenase 8) (EC 1.1.1.239)

Note

For Research Use Only , Not for Diagnostic Use