


Product Name: Recombinant Human ELMO1 Protein
Catalog #: PEH7837



Summary

Name	ELMO1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q92556
Host	E.coli
Species	Human
Molecular Mass	28 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	QFKSKLQNLS YTEILKIRQS ERMNQEDFQS RPILELKEKI QPEILELIKQ QRLNRLVEGT CFRKLNARRR QDKFWYCRLS PNHKVLHYGD LEESPQGEVP HDSLQDKLPV ADIKAVVTGK DCPHMKEKGA LKQNKEVLEL AFSILYDSNC QLNFIAPDKH EYCIWTDGLN ALLGKDMMSD LTRNDLDTLL SMEIKLRLLD LENIQIPDAP PPIPKESNY

image



Background

Alternative Names Engulfment and cell motility protein 1 (Protein ced-12 homolog)

Note

Product Name: Recombinant Human ELMO1 Protein
Catalog #: PEH7837



For Research Use Only , Not for Diagnostic Use