

Product Name: Recombinant Human NR4A3 Protein
Catalog #: PEH7838



Summary

Name	NR4A3_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q92570
Host	E.coli
Species	Human
Molecular Mass	11 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	TCAVCGDNAA CQHYGVRTCE GCKGFFKRTV QKNAKYVCLA NKNCPVDKRR RNRCQYCRFQ KCLSVGMVKE VVRTDSLKGR

image



Background

Alternative Names	Nuclear receptor subfamily 4 group A member 3 (Mitogen-induced nuclear orphan receptor) (Neuron-derived orphan receptor 1) (Nuclear hormone receptor NOR-1)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use