


Product Name: Recombinant Human CREBBP Protein
Catalog #: PEH7842



Summary

Name	CBP_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-MBP & C-His
Accession #	Q92793
Host	E.coli
Species	Human
Molecular Mass	74 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	PTPGSVPSAT QTQSTPTVQA AAQAQVTPQP QTPVQPPSVA TPQSSQQQPT PVHAQPPGTP LSQAAASIDN RVPTPSSVAS AETNSQQPGP DVPVLEMKTE TQAEDTEPDP GESKGEPRSE MMEEDLQGAS QVKEETDIAE QKSEPMEVDE KKPEVKVEVK EEESSSNGT ASQSTSPSQP RKKIFKPEEL RQALMPTLEA LYRQDPESLP FRQPVDPQLL GIPDYFDIVK NPMDLSTIKR KLDTGQYQEP WQYVDDVWLM FNNAWLYNRK

image



Background

Alternative Names	CREB-binding protein (Histone lysine acetyltransferase CREBBP) (EC 2.3.1.48) (Protein-lysine acetyltransferase CREBBP) (EC 2.3.1.-)
--------------------------	--

Product Name: Recombinant Human CREBBP Protein
Catalog #: PEH7842



Note

For Research Use Only , Not for Diagnostic Use