


**Product Name: Recombinant Human INPP5D Protein**  
**Catalog #: PEH7844**



## Summary

<b>Name</b>	SHIP1_HUMAN
<b>Purity</b>	Greater than 85% as determined by reducing SDS-PAGE
<b>Endotoxin level</b>	Please contact with the lab for this information
<b>Tag</b>	C-His
<b>Accession #</b>	Q92835
<b>Host</b>	E.coli
<b>Species</b>	Human
<b>Molecular Mass</b>	13 kDa
<b>Buffer</b>	Phosphate buffered saline
<b>Shipping</b>	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
<b>Stability&amp;Storage</b>	Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
<b>Sequence</b>	VKSQEGENEE GSEGELVVKF GETLPKPKPI ISDPEYLLDQ HILISIKSSD SDESYGEGCI ALRLEATETQ LPIYTPLTHH GELTGHFQGE IKLQTSQGKT

## image



## Background

<b>Alternative Names</b>	Phosphatidylinositol 3,4,5-trisphosphate 5-phosphatase 1 (EC 3.1.3.86) (Inositol polyphosphate-5-phosphatase D) (EC 3.1.3.56) (Inositol polyphosphate-5-phosphatase of 145 kDa) (SIP-145) (Phosphatidylinositol 4,5-bisphosphate 5-phosphatase) (EC 3.1.3.36) (SH2 domain-containing inositol 5'-phosphatase 1) (SH2 domain-containing inositol phosphatase 1) (SHIP-1) (p150Ship) (hp51CN)
--------------------------	---

**Product Name: Recombinant Human INPP5D Protein**  
**Catalog #: PEH7844**



---

**Note**

For Research Use Only , Not for Diagnostic Use