


Product Name: Recombinant Human ARHGEF1 Protein
Catalog #: PEH7845



Summary

Name	ARHG1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q92888
Host	E.coli
Species	Human
Molecular Mass	27 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	PGPSRPGLVP VSIIGAEDD FENELETNSE EQNSQFQSLE QVKRRPAHLM ALLQHVALQF EPGPLLCLH ADMLGSLGPK EAKKAFLDFY HSFLEKTAVL RVPVPPNVAF ELDRTRADLI SEDVQRRFVQ EVVQSQQVAV GRQLEDFRSK RLMGMTPEEQ ELAQLEAWVG RDRASYEARE RHVAERLLMH LEEMQHTIST DEEKSAVVN AIGLYMRHLG

image



Background

Alternative Names	Rho guanine nucleotide exchange factor 1 (115 kDa guanine nucleotide exchange factor) (p115-RhoGEF) (p115RhoGEF) (Sub1.5)
--------------------------	---

Note

Product Name: Recombinant Human ARHGEF1 Protein
Catalog #: PEH7845



For Research Use Only , Not for Diagnostic Use