


Product Name: Recombinant Human CUL5 Protein
Catalog #: PEH7848



Summary

Name	CUL5_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q93034
Host	E.coli
Species	Human
Molecular Mass	40 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	FDLFSDVHAV CLWDDKGPAP IHQALKEDIL EFIKQAQARV LSHQDDTALL KAYIVEWRKF FTQCDILPKP FCQLEITLMG KQGSNKKSNV EDSIVRKLML DTWNESIFSN IKNRLQDSAM KLVHAERLGE AFDSQLVIGV RESYVNLCSN PEDKQLIYRD NFEKAYLDST ERFYRTQAPS YLQQNGVQNY MKYADAKLKE EEKRALRYLE TRRECNSVEA LMECCVNALV TSFKETILAE CQGMIKRNET EKLHLMFSLM DKVPNGIEPM LKDL EEHIS AGLADMVAAA ETITDSEKY VEQLLTLFNR FSKLVKEAFQ DDPFRFLTARD

image



Background

Alternative Names Cullin-5 (CUL-5) (Vasopressin-activated calcium-mobilizing receptor 1) (VACM-1)

Product Name: Recombinant Human CUL5 Protein
Catalog #: PEH7848



Note

For Research Use Only , Not for Diagnostic Use