



**Product Name: Recombinant Human PBXIP1 Protein**  
**Catalog #: PEH7850**



## Summary

<b>Name</b>	PBIP1_HUMAN
<b>Purity</b>	Greater than 85% as determined by reducing SDS-PAGE
<b>Endotoxin level</b>	Please contact with the lab for this information
<b>Tag</b>	C-His
<b>Accession #</b>	Q96AQ6
<b>Host</b>	E.coli
<b>Species</b>	Human
<b>Molecular Mass</b>	45 kDa
<b>Buffer</b>	Phosphate buffered saline
<b>Shipping</b>	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
<b>Stability&amp;Storage</b>	Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
<b>Sequence</b>	GVCLSGGRGP QGDKAIREQG PREQEPELSF LKQKEQLEAE AQALRQELER QRRLLGSVQQ DLERSLQDAS RGDPAHAGLA ELGHRLAQLK QGLENWGQDP GVSANASKAW HQKSHFQNSR EWSGKEKWWD GQRDRKAEHW KHKKEESGRE RKKNWGGQED REPAGRWKEG RPRVEESGSK KEGKRQGPKE PPRKSGSFHS SGEKQKQPRW REGTKDSDHP LPSWAELLRP KYRAPQGCSG VDECARQEGL TFFGTELAPV RQQELASLLR TYLARLPWAG QLTKELPLSP AFFGEDGIFR HDRLRFRDFV DALEDSLEEV AVQQTGDDDE VDDFEDFIFS HFFGDKALKK RSGKKDKHSQ SPRAAGPREG HSHSHHHHHR

## image



## Background

**Alternative Names** Pre-B-cell leukemia transcription factor-interacting protein 1 (Hematopoietic

**Product Name: Recombinant Human PBXIP1 Protein**  
**Catalog #: PEH7850**



---

PBX-interacting protein)

**Note**

For Research Use Only , Not for Diagnostic Use