



Product Name: Recombinant Human DHX58 Protein
Catalog #: PEH7853



Summary

Name	DHX58_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q96C10
Host	E.coli
Species	Human
Molecular Mass	24 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	FDDRKNELAH LATHGPENPK LEMLEKILQR QFSSSNSPRG IIFTRTRQSA HSLLLWLQQQ QGLQTVDIRA QLLIGAGNSS QSTHMTQRDQ QEVIQKFQDG TLNLLVATSV AEGLDIPHC NVVRYGLLT NEISMVQARG RARADQSVYA FVATEGSREL KRELINEALE TLMEQAVAAV QKMDQAEYQA KIRDLQQAAL

image



Background

Alternative Names	Probable ATP-dependent RNA helicase DHX58 (EC 3.6.4.13) (Probable ATP-dependent helicase LGP2) (Protein D11Lgp2 homolog) (RIG-I-like receptor 3) (RLR-3) (RIG-I-like receptor LGP2) (RLR)
--------------------------	---

Note

Product Name: Recombinant Human DHX58 Protein
Catalog #: PEH7853



For Research Use Only , Not for Diagnostic Use