

Product Name: Recombinant Human HOOK2 Protein
Catalog #: PEH7857



Summary

Name	HOOK2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q96ED9
Host	E.coli
Species	Human
Molecular Mass	39 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	HIQRIMTLEE SVQHVVMEAI QELMTKDTPD SLSPETYGNF DSQSRRYYFL SEEAEEGDEL QQRCLDLERQ LMLLSEEKQS LAQENAGLRE RMGRPEGEGT PGLTAKKLLL LQSQLEQLQE ENFRLESGRE DERLRCAELE REVAELQHRN QALTSLAQEA QALKDEMDEL RQSSERAGQL EATLTSCRRR LGELRELRRQ VRQLEERNAG HAERTRQLED ELRRAGSLRA QLEAQRQVQ ELQGQRQEEA MKAEKWLFEC RNLEEKYESV TKEKERLLAE RDSLREANEE LRCAQLQPRG LTQADPSLDP TSTPVDNLAA

image



Background

Alternative Names Protein Hook homolog 2 (h-hook2) (hHK2)

Product Name: Recombinant Human HOOK2 Protein
Catalog #: PEH7857



Note

For Research Use Only , Not for Diagnostic Use