


Product Name: Recombinant Human MOCOS Protein
Catalog #: PEH7858



Summary

Name	MOCOS_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q96EN8
Host	E.coli
Species	Human
Molecular Mass	22 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	SQRNAKKKHG KDQLPGTMAT LSLVNEAQL LINTSSILEL HRQLNTSDEN GKEELFSLKD LSLRFRANII INGKRAFEEE KWDEISIGSL RFQVLGPCHR CQMICIDQQT GQRNQHVFK LSEPRETKVN FGMVLMHASL DLSSPCFLSV GSQVLPVLKE NVEGHDLPAS

image



Background

Alternative Names Molybdenum cofactor sulfurase (MCS) (MOS) (MoCo sulfurase) (hMCS) (EC 2.8.1.9) (Molybdenum cofactor sulfurtransferase)

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human MOCOS Protein
Catalog #: PEH7858

