


Product Name: Recombinant Human PGM2 Protein
Catalog #: PEH7859



Summary

Name	PGM2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q96G03
Host	E.coli
Species	Human
Molecular Mass	38 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	EDARLDQETA QWLRWDKNSL TLEAVKRLIA EGNKEELRKC FGARMEFGTA GLRAAMGPGI SRMNDLTIQ TTQGFCRYLE KQFSDLKQKG IVISFDARAH PSSGGSSRRF ARLAATTFIS QGIPVYLFSD ITPTFPVPFT VSHLKLKAGI MITASHNPKQ DNGYKVYWDN GAQIISPHDK GISQAIEENL EPWPQAWDDS LIDSSPLLHN PSASINNDYF EDLKKYCFHR SVNRETKVKF VHTSVHGVGH SFVQSAFKAF DLVPPEAVPE QKDPDPEFPT VKYPNPEEGK GVLTLSFALA DKTKARIVLA NDPDADRLAV

image



Background

Alternative Names	Phosphoglucomutase-2 (PGM 2) (EC 5.4.2.2) (Glucose phosphomutase 2) (Phosphodeoxyribomutase) (Phosphopentomutase) (EC 5.4.2.7)
--------------------------	--

Product Name: Recombinant Human PGM2 Protein
Catalog #: PEH7859



Note

For Research Use Only , Not for Diagnostic Use