


**Product Name: Recombinant Human PDLIM5 Protein**  
**Catalog #: PEH7861**



## Summary

<b>Name</b>	PDLI5_HUMAN
<b>Purity</b>	Greater than 85% as determined by reducing SDS-PAGE
<b>Endotoxin level</b>	Please contact with the lab for this information
<b>Tag</b>	C-His
<b>Accession #</b>	Q96HC4
<b>Host</b>	E.coli
<b>Species</b>	Human
<b>Molecular Mass</b>	37 kDa
<b>Buffer</b>	Phosphate buffered saline
<b>Shipping</b>	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
<b>Stability&amp;Storage</b>	Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
<b>Sequence</b>	PAPWGFRLQG GKDFNMPLTI SSLKDGKAA QANVRIGDVV LSIDGINAQQ MTHLEAQNKI KGCTGSLNMT LQRASAAPKP EPVPVQKGEP KEVVKVPVIT SPAVSKVTST NNMAYNKAPR PFGSVSSPKV TSIPSPSSAF TPAHATTSSH ASPSPVAAVT PPLFAASGLH ANANLSADQS PSALSAGKTA VNVPRQPTVT SVCSETSQEL AEGQRRGSQG DSKQQNGPPR KHIVERYTEF YHVPTHSDAS KKRLIEDTED WRPRTGTTQS RSFRILAQIT GTEHLKESEA DNTKKANNSQ EPSPQLASSV ASTRSMPEL DSPTSGRPGV

## image



## Background

<b>Alternative Names</b>	PDZ and LIM domain protein 5 (Enigma homolog) (Enigma-like PDZ and LIM domains protein)
--------------------------	---

**Product Name: Recombinant Human PDLIM5 Protein**  
**Catalog #: PEH7861**

---



**Note**

For Research Use Only , Not for Diagnostic Use