


Product Name: Recombinant Human ITCH Protein
Catalog #: PEH7862



Summary

Name	ITCH_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q96J02
Host	E.coli
Species	Human
Molecular Mass	22 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	KDETRVSTNG SDDPEDAGAG ENRRVSGNNS PLSNNGGFKP SRPPRPSRPP PPTPRRPASV NGSPSATSES DGSSTGSLPP TNTNTNTSEG ATSGLIPLT ISGGSGPRPL NPVTQAPLPP GWEQRVDQHG RYYYVDHVEK RTTWDRPEPL PPGWERRVDN MGRIYYVDHF TRTTTWQRPT

image



Background

Alternative Names	E3 ubiquitin-protein ligase Itchy homolog (Itch) (EC 2.3.2.26) (Atrophin-1-interacting protein 4) (AIP4) (HECT-type E3 ubiquitin transferase Itchy homolog) (NFE2-associated polypeptide 1) (NAPP1)
--------------------------	---

Note

Product Name: Recombinant Human ITCH Protein
Catalog #: PEH7862



For Research Use Only , Not for Diagnostic Use