


Summary

Name	HUTU_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q96N76
Host	E.coli
Species	Human
Molecular Mass	33 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	DSLDRCIQRL REARKKKEVL SLGYHGNVVA LWERLVHELD TTGECLVDLG SDQTSCHNPF NGGYYPVQLS FTEAQLMAS NPAVFKDLVQ ESLRRQVSAI NRLAEEKFFF WDYGNAFLLE AQRAGADVEK KGAGRTEFRY PSYVQHIMGD IFSQGFDPFR WVCTSGDPQD LAVTDELATS VLEEAIADGV KSVKQLQYMD NIRWIREAAR HRLVVGSQLR ILYSDQKGRV AIAVAINQAI ACRRIKAPVV LSRDHHDVSG TDSPFRETSN

image



Background

Alternative Names	Urocanate hydratase (Urocanase) (EC 4.2.1.49) (Imidazolonepropionate hydrolase)
--------------------------	---

Product Name: Recombinant Human UROC1 Protein
Catalog #: PEH7865



Note

For Research Use Only , Not for Diagnostic Use