


**Product Name: Recombinant Human GBP5 Protein**  
**Catalog #: PEH7868**



## Summary

<b>Name</b>	GBP5_HUMAN
<b>Purity</b>	Greater than 85% as determined by reducing SDS-PAGE
<b>Endotoxin level</b>	Please contact with the lab for this information
<b>Tag</b>	C-His
<b>Accession #</b>	Q96PP8
<b>Host</b>	E.coli
<b>Species</b>	Human
<b>Molecular Mass</b>	29 kDa
<b>Buffer</b>	Phosphate buffered saline
<b>Shipping</b>	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
<b>Stability&amp;Storage</b>	Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
<b>Sequence</b>	VQLPMETLQE LLDLHRTSER EAIEVFMKNS FKDVDQSFQK ELETLLDAKQ NDICKRNLEA SSDYCSALLK DIFGPLEEAV KQGIYSKPGG HNLFIQKTEE LKAKYYREPR KGIQAEVLQ KYLKSSESVS HAILQTDQAL TETEKKKKEA QVKAEAEKAE AQRLAAIQRQ NEQMMQERER LHQEQVRQME IAKQNWLAEQ QKMQEQQMQE QAAQLSTTFQ AQNRSLLSEL QHAQRTVNND

## image



## Background

**Alternative Names** Guanylate-binding protein 5 (EC 3.6.5.-) (GBP-TA antigen) (GTP-binding protein 5) (GBP-5) (Guanine nucleotide-binding protein 5)

## Note

**Product Name: Recombinant Human GBP5 Protein**  
**Catalog #: PEH7868**



---

For Research Use Only , Not for Diagnostic Use